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
P: 800-321-0607
F: 800-329-2734

AW-1.02070 REV J 0311


CAUTION: FEDERAL LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A LICENSED HEALTH CARE PRACTITIONER.
CAUTION: FOR SINGLE PATIENT USE ONLY.
CAUTION: IF YOU EXPERIENCE INCREASED PAIN, SWELLING, OR ANY ADVERSE REACTIONS WHILE USING THIS PRODUCT, IMMEDIATELY CONSULT YOUR MEDICAL PROFESSIONAL.
WARNING: CAREFULLY READ USE INSTRUCTIONS AND WARNINGS PRIOR TO OPERATION.
WARNING: THIS DEVICE WILL NOT PREVENT OR REDUCE ALL INJURIES. PROPER REHABILITATION AND ACTIVITY MODIFICATION ARE ALSO AN ESSENTIAL PART OF A SAFE TREATMENT PROGRAM

ROADRUNNER KNEE BRACES


ROADRUNNER

SHORTRUNNER

ROADRUNNER
Wrap-On


SHORTRUNNER
Wrap-On


APPLICATION INSTRUCTIONS FOR KNEE BRACES:

The wrap-on brace is applied with the closures in the front of the leg and the kneecap cutout centered over the knee. Secure the closures firmly.



The slip-on brace is applied by pulling the brace on with the kneecap cutout centered over the knee.

Both the wrap-on and pull-on braces have tightening straps. The straps should be tightened firmly, starting with the two straps closest to the knee.

ADJUSTMENT INSTRUCTIONS
Polycentric hinge with stops

Both flexion and extension stops are supplied with the brace. The brace comes with a 10° extension stop already installed. The stops are labeled in degrees as follows;

Extension Stops: 0°, 10°, 20°, 30° & 40°
Flexion Stops: 45°, 60°, 75° & 90°
(no flexion stops are installed upon shipment)



To change the stops, first, completely remove the screw as shown. Then remove the stop. Insert the new stop, align with the holes in the hinge, then reinsert the screw. Using a Phillips head screw driver, ensure that the screws are seated in the recess and completely tightened. Both medial (inside) and lateral (outside) hinge stops must be identical in degree.

NOTE: Due to soft tissue compression, we recommend the use of extension stops at least 10° greater than the desired limit. (e.g., if you desire a 10° extension limit, use the 20° extension stop).

TIPS ON USE AND CARE OF YOUR KNEE BRACE:

- After exposure to salt water or dirt, thoroughly rinse the brace with fresh water and allow to air dry.
- A dry lubricant, such as teflon spray, may be used on the hinge occasionally. Do not use wet lubricants.
- The hinge assemblies may be removed to allow easier washing of the sleeve and straps.
- The side bars may be bent to accommodate individual legs. Bend in equal amounts above and below the hinge.
- Periodically inspect the sleeve, straps and hinges for excessive wear. Be sure all screws are tight and stops in place.

ADJUSTMENT INSTRUCTIONS
Polycentric hinge with pins

The hinge settings are labeled in degrees as follows;

Extension settings: 0°, 10°, 30°, 50°, 70°, 90° & 110°
Flexion settings: 0°, 10°, 30°, 50°, 70°, 90° & 110°

To change the settings, first, open the hinge cover as shown. Then, move the pins to the desired extension and/or flexion settings for both the upper and lower leg. When satisfied, close the hinge cover. Both medial (inside) and lateral (outside) hinge settings must be identical in degree.

NOTE: Due to soft tissue compression, we recommend the use of extension stops at least 10° greater than the desired limit. (e.g., if you desire a 10° extension limit, use the 20° extension stop).

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